

A High Efficiency Instrument for Collecting Airborne Particles Down to 10 Nanometers on Solid Surfaces or Into Liquid

Presented at Healthy Buildings America, Boulder, CO, July 19-22, 2015

Arantzazu EIGUREN-FERNANDEZ^{1*}, Nick BUGLAK², Ning LI², Pat KEADY³, and Susanne HERING^{1,3}
¹Aerosol Dynamics Inc., Berkeley, CA; ² Dept. of Pathobiology & Diagnostic Investigation, Michigan State Univ., East Lansing, MI; ³Aerosol Devices Inc, Fort Collins, CO



INTRODUCTION

Airborne microorganisms (viruses, bacteria, and fungi) are part of ambient particulate matter found as either individual particles, aggregates, or bound to other airborne particles. Their size distribution ranges from a few nanometers to several micrometers.

Active sampling methods with liquid collectors (best for maintaining microbe viability), filters, and impactors¹ are widely used. However, they are limited by: *inefficient physical collection*, <10% for particles 30-100nm² and *low viability of those collected*, due to microbe stress, desiccation and breakup.

We present an alternative method of sampling airborne particles that uses water condensational growth and soft impaction on either a dry surface or in liquid media. Our **SPOT SAMPLER™** collector provides

high physical collection efficiencies for particles 10 nm and larger, and **higher viability** of the collected organisms

INSTRUMENT CAPABILITIES

The Spot Sampler system collects particles via gentle impaction on a solid substrate or as suspensions into liquid for chemical, biological and toxicological characterization^{3,4,5} (Figure 1).




Figure 1. a) Integrated Spot Sampler; b) Liquid collection; c) Dry multi-well collection disk

Advantages of aerosol particle collection into liquids (Figure 1b):

- **Direct collection** into liquid reduces sampling artifacts and maintains microorganism viability
- **No need** for extraction and resuspension steps
- Collects **soluble particle-phase** components, and **insoluble particles** suspended in the collection liquid
- The **concentrated nature** of the collection allows reducing sampling times
- **Versatile collection medium**: water, culture media, virus growth medium, etc.
- Provides concentrated particle samples ready to be characterized using **rapid microbiological analysis** techniques, such as PCR, and **in-vitro assays**⁵

Advantages of dry particle collection (Figure 1a, 1c):

- **Uninterrupted, time-resolved** collection of concentrated “spots” (1-mm) in a 33-well disk: minutes to hours
- **Easy transport** of sample plate from the field to the laboratory for analysis
- **Automated** extraction and injection by autosampler reduces contamination for wet chemical analysis^{3,4} (IC, HPLC, etc.)
- **Concentrated deposits** are suitable for spectrographic or molecular analysis

HOW IT WORKS

The sampler uses a three-stage, moderated laminar-flow condensation method to grow tiny airborne particles as small as 5nm into ~3µm droplets and collects them by bounce-free, soft impaction into the collection substrate^{3,6} (Figure 2).




Figure 2. Three-stage water condensational growth approach




Figure 3. Temperature, saturation ratio, and droplet size under standard operating conditions⁶

The moderated temperatures under 30°C (Figure 3), create a hospitable environment advantageous for viable microorganism collection. The time the particle is subjected to high humidity conditions is brief – on the order of milliseconds. If the presence of water is an issue for the particles or assay, a dry deposit eliminates any negative effects due to condensation.

PARTICLE COLLECTION EFFICIENCY




Figure 4 Collection efficiency test with monodisperse hydrophilic and hydrophobic aerosols

The physical collection efficiency is over 90 percent for liquid and dry collection down to 5nm for hydrophilic particles and 10nm for highly hydrophobic particles.

IN-VITRO TOXICOLOGY ANALYSIS

In-vitro toxicity of airborne particles (promising early results):

- Collect ambient particles as a liquid suspension over a period of four days at Michigan State University (East Lansing, MI)
- Expose human bronchial epithelial **cells line BEAS-2B** to the particle suspension for 16 hours
- Measure pro-inflammatory markers interleukin-6 (IL-6) and interleukin-8 (IL-8)
- Ambient aerosols induced a 10- to 34-fold increase in IL-6, and a 5- to 25-fold increase for IL-8, compared to controls (Figure 5)




Figure 5. Production of IL-6 and IL-8 in cells after exposure to East Lansing ambient aerosols.

CONCLUSIONS

1. The Spot Sampler instrument provides a new approach to efficiently collect airborne particles as concentrated, ready-to-analyze, liquid suspensions or dry deposits with low contamination.
2. High physical collection efficiencies (>90% from 10nm to 10µm) for both collection configurations, increasing the range of efficient collection from nm (bare viruses) to µm (virus associated with airborne particles, bacteria and fungi).
3. Gentle collection, moderate temperatures, and the ability to control the sampling relative humidity conditions, is advantageous for collecting viable microorganisms.

REFERENCES

1. Verrault D., Moineau S., Duchaine C. - Methods for Sampling of Airborne Viruses. Microbiology and Molecular Biology Reviews (2008) 72(2): 413-444
2. Hogan C.J., Jr., E.M. Kettleson, M.H. Lee, B. Ramaswami, L.T. Angenent, and P. Biswas - Sampling Methodologies and dosage assessment techniques for submicrometre and ultrafine virus aerosol particles. *J. Appl. Microbiol.* (2005) 99:1422-1434
3. Eiguren Fernandez, A., Lewis G.S., and Hering, S.V – Design and Laboratory Evaluation of a Spot Sampler for Time-Resolved Measurement of Airborne Particle Composition. *Aerosol Sci. Tech.* (2014) 48(6):655-663
4. Eiguren Fernandez A., Lewis G.S., Spielman S et al. - Time-resolved characterization of particle associated polycyclic aromatic hydrocarbons using a newly-developed sequential spot sampler with automated extraction and analysis. *Atmos. Environm.* (2014) 96:125-134
5. Lewis G., Buglak N., Spielman S. et al - Concentrated Aerosol Collection into Liquid for Toxicological Studies. American Assoc. for Aerosol Research, Orlando, FL, 20-24 Oct. 2014
6. Hering S, Lewis G, Spielman S. - Moderated, Water-Based Condensational Particle Growth in a Laminar Flow. *Aerosol Sci. Tech.* (2014) 48(4):401-8

ACKNOWLEDGEMENTS

This study was made possible by the support from National Institutes of Health (Grants: RC3 ES019081-01 & R43-ES022523). Authors would like to thank Dr. Ryan Lewandowski and Paula Palmiter from Michigan State University for their help during field deployment. Collector technology is licensed by Aerosol Devices Inc. from Aerosol Dynamics Inc. with U.S. Patents #6712881, #7736421, #8801838 and German Patent #10392241. Other patents pending.